

TECHNICAL SPECIFICATIONS : Electrical muscle stimulator

1. It has advantage of manual selection combined with the technology of embedded microcontroller resulting in efficient usage
2. Timer: Can be used to set the treatment time to a maximum of 30 minutes
3. Intensity control can be used to select continuously 0.01 to 50.0 ma and controlled with manual controller separately.
4. Displays the mode of treatment as F and G for Faradic and Galvanic respectively while using the program selector.
5. Output terminals: Red and Black Terminals for connecting the patient or output cable.
6. Red round led between the output terminals will illuminate corresponding to the duration of pulses of both F7G outputs